

HF440

CRRT

The All in One Device



HF₄₄₀ offers the most complete range of treatments

- **SCUF:** Slow Continuous UltraFiltration
 - **CVVH:** Continuous Veno-Venous Hemofiltration
 - **CVVHD:** Continuous VenoVenous HemoDialysis
 - **CVVHDF:** Continuous VenoVenous HemoDiaFiltration
 - **HVHF:** Continuous VenoVenous HemoDiaFiltration
 - **TPE:** Therapeutic Plasma Exchange
 - **DFPP:** Cascade filtration Double Filtration
 - **HP:** Hemoperfusion
 - **PP:** Plasmaperfusion
 - **CPFA:** Coupled Plasma Filtration Adsorption
-
- HF440 offers the most complete range of treatments, all available with a full range of flows from small babies to large adults. Single needle blood access, as well as all anticoagulation methods are also available for all treatments, thus the HF440 perfectly meets all the hospital needs.
 - The HF440 is easy to use thanks to its touchscreen interface, the record of data and the graphs for pressures and weights. Its low weight and small dimensions make it easy to move while the pre-assemble sets, scales system up to 40 kg and clear messages helps reducing operators workload.
 - HF440 offers unique features such as automatic adjustment of predilution, rinsing by flush of the cascade filter, single needle blood access or CPFA treatment.
 - Baby to adult use with adapted flows, volumes and safety


HF440

Full touch screen

The HF440 is easy to use thanks to its full touch screen interface, the record of data and the graphs for pressures and masses

Though for operators

Its light weight and small dimensions make it easy to move while the pre-assembled sets, scales system up to 40 kg and clear messages helps to reduce operators workload



Scales For up to 8 bags and 40 liters for substitution, dialysate and waste.

Heater Integrated plate fluid heater, temperature setting 35-38°C, temperature graph overtime

Technical overview

Touchscreen, 5 pumps, 3 scales up to 8 bags of 5L each, heater, syringe pump also available during rinsing, save rinsing, automatic predilution, pre and postdilution ratio editable during the treatment or automatically set, air detector, blood leak detector, single needle activable during treatment, sound and visual alarm, on screen alarm guidance, profiles saving, more than 20 minutes backup battery, care time button, wi-fi, bluetooth and more.

Please contact your local representative for more information

Manufacturing Plant - Factory:

Nipro India Corporation Pvt. Ltd. Plot No. E-1/1, Khandala Industrial Area PH-1, Taluka - Khandala, District-Satara, Maharashtra, India- Pin Code- 412 802. Tel.: +91 2169 3050 50/ 3072 00 Email: indiaho@nipromedindia.com

Sales Head Office:

Nipro Medical India Pvt. Ltd. GMR Hyderabad Int'l Airport, 2nd floor, GMR Aero Tower, Shamshabad, Hyderabad - 500409. Andhra Pradesh. Tel.: +91 40 2459 9499. Email: salesho@nipromedindia.com

Kolkata Office:

Nipro Medical India Pvt. Ltd. Office number 107, First floor, C Block, DLF IT park, 08 - Major Arterial Road, Block - AF, New Town Rajarhat, Kolkata - 700 156, West Bengal. Tel.: +91 33 4067 9115 / 19 Telefax: +91 33 4067 9114. Email: eastbo@nipromedindia.com

Mumbai Office:

Nipro Medical India Pvt. Ltd. E-10, Unit No. 1, 2 & 3, Bhumi World Industrial Park, Pimpas Village, Mumbai-Nashik Highway, Before Kalyan-Bhiwandi Naka, Opp Tata Amantra, Pimpas, Thane, Maharashtra 421302. Tel.: +91 22 66114545. Email: westcom@nipromedindia.com

Gurgaon Office:

Nipro Medical India Pvt. Ltd. Unit No. 306, 3rd Floor, Time Tower, M G Road, Gurgaon - 122 002, Haryana. Tel.: +91 124 427 8601/437 8601. Telefax: +91 124 437 8601. Email: northbo@nipromedindia.com

Chennai Office:

Nipro Medical India Pvt. Ltd. Flat no. 59, 4th floor, Krishna Building, Building no. -41, 1st avenue, 100 feet road, Near Hot chips restaurant and Lakshmi Mahal, Ashok Nagar, Chennai - 600 089. Tamilnadu. Tel.: +91 44 24895560 / 70. Mobile: +91 96291 7682. Email: southbo@nipromedindia.com

Ahmedabad Office:

Nipro Medical India Pvt. Ltd. Sukh Shanti complex, Ground floor Hall, Opp. Jekson Hydraulics, Near Bhagyodaya Hotel, Chandodhar, Taluka-Sanand, Ahmedabad - 382213, Gujarat. Tel. No. : +91 9426282189. Email: ahmedabadbo@nipromedindia.com

Lucknow Office:

Nipro Medical India Pvt. Ltd. E-25, Phase -11, Transport Nagar, Lucknow -226012, Uttar Pradesh. Tel. No. : +91 7860537725. Email: lucknowbo@nipromedindia.com

Jaipur Office:

Nipro Medical India Pvt. Ltd. F-1291, SITAPURA PH-III, SITAPURA INDUSTRIAL AREA, JAIPUR-302022, Rajasthan. Tel. No. : +91 141271769. Email: jaipurbo@nipromedindia.com

Mohali Office:

Nipro Medical India Pvt. Ltd. 703, GLPL, Sector No.-82, Mohali - 160055, Punjab. Tel. No. : Email: mohalibo@nipromedindia.com

Indore Office:

Nipro Medical India Pvt. Ltd. Shop No/Plot No. B-16, G.F. And First Floor, Sch. No. 78, Part 1, New Loha Mandi, Tehsil: Indore, Dist: Indore, Pin: 452010, Madhya Pradesh. Tel. No. : +91 9407051912. Email: indorebo@nipromedindia.com

Raipur Office:

Nipro Medical India Pvt. Ltd. Shop. No. 73 & 74, Aushadi Vatika, Medical Complex Dumartari, New Dhantari Road, Block-Dharsinwa Raipur - 492015, Chhattisgarh. Tel. No. : +91 7716900070. Email: raipurbo@nipromedindia.com

Bangalore Office:

Nipro Medical India Pvt. Ltd. SY No.51/1 A2, Hegadadevanapura, Opp. Unibic Factory, Huskur Road, Dasanapura, Hubli, Bangalore - 562123, Karnataka. Tel. No. : +91 80 23718359/60. Email: bangalorebo@nipromedindia.com

Cochin Office:

Nipro Medical India Pvt. Ltd. 16/1D, Anjith Arcade, Kaippadamugal, Seaport Airport Road, Vadakode Po, Kallamssery, Cochin - 682021, Kerala. Tel. No. : +91 7736062123. Email: cochinbo@nipromedindia.com

NIPRO DIALYSIS TREATMENT & TRAINING CENTER Office

No.301, 306 & 307, 3rd Floor, Lalani Quantum, Above Chellaram Multi-specialty Hospital, Mumbai-Pune Bypass Rd, Patil Nagar, Bavdhan, Pune - 411 021 Maharashtra. Tel.: +91 20-67928700. Email: punebo@nipromedindia.com. To book appointment for Dialysis Treatment and Dialysis Training Course pls contact 020-67928700.